

How did we measure a floating "2 µm cell" on a chip?

Synechocystis sp. PCC6803

 Floating unicellular cyanobacterium Model microorganism: photosynthesis, biofuel and acclimation to environmental changes.









Wild type

Mutant type ($\Delta MscL$)

Purpose : Evaluate the mechanical properties of wild type and mutant type cell under osmotic shock.

Direct drive of on-chip probe, and Sampling moiré method





Contact person : Shinya SAKUMA E-mail: sakuma@mech.nagoya-u.ac.jp URL: http://www.biorobotics.mech.nagoya-u.ac.jp/ Dept. of Micro/nano System Engineering, Nagoya University TEL: +81-52-789-5220, FAX : +81-52-789-5027

We report the evaluation results of mechanical characterization of Synechocystis sp. PCC 6803. We constructed the force measurement system using the microfluidic chip. Using the system, we evaluated the Young's modulus of two cell groups of wild type and mutant type, which was knocked out mechanosensitive channels. We confirmed that We succeeded in measuring the mechanical characteristics of cyanobacterium.

